

Technical details and specification

Computer unit PC: Compaq PC with 15" TFT screen ∷ Notebook: Compaq notebook

DAQ board PC: 16 channels: 200kS/s, 12 Bit, 8 digital I/Os, 0 or 2 analog outputs
 16 channels: 100kS/s, 16 Bit, 8 digital I/Os, 0 or 2 analog outputs
 64 channels: 500kS/s, 12 Bit, 8 digital I/Os, 0 or 2 analog outputs

 Notebook: 16 channels: 200kS/s, 12 Bit, 8 digital I/Os, 2 analog outputs
 16 channels: 20kS/s, 16 Bit, 8 digital I/Os

Connector box input/output front side: 2 pin Cinch connector to analog inputs, digital inputs/outputs, analog outputs, +5v
 rear side: 25 pin Sub-D connector to *g.BSamp*, 9pin Sub-D connector to digital inputs/outputs, analog outputs, +5v
 User-specific connector boxes upon request

Dimension 16 channel version: length: 31 cm, width: 25 cm, height: 7 cm
 32-64 channel version: length: 31 cm, width: 25 cm, height: 10 cm